Balloon-Lift-Control

developed with Prof. Dr. Benner, Dr. Dr. Heuckmann, Dr. Bauer

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
- Minimally invasive lifting of internal sinus
- Optimal use of existing bone substance independent of vertical height
- Safe, defined and simple access to the maxillary sinus mucosa
- Complication free lifting by up to approx. 10 mm through extended balloon
- The procedure is easy to pick up and handle
- The use of the procedure saves time and costs
- Suitable to prepare the jaw bone for the implantation of all current systems
- Possibility of the direct implantation by the application of the Osteotoms I with the diameter 3.80 mm

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
¹ 204=RA, 205=RA L, 206=RA XL

² Largest working part diameter in 1/10 mm

Minimaldurchmesser Minimal diameter Diámetro mínimo

Außendurchmesser External diameter Diámetro externo

Führungshülsen für Osteotom I
Guidance shells for osteotome I
Manguitos de guía de osteótomo I

CE 0044

Fig.	BLC0B	BLC0B	BLC0B	BLC0B	BLC0B
Größe Size ² Tamaño	045	045	045	045	045
Länge Length	9,0	10,0	11,0	12,0	13,0
Longitud mm					
Bohrtiefe Drilling depth	5,0	4,0	3,0	2,0	1,0
Profundidad de perforación mm					

¹ 204=RA, 205=RA L, 206=RA XL

² Largest working part diameter in 1/10 mm

Führungshülsen für Osteotom II
Guidance shells for osteotome II
Maguito de guía de osteótomo II



Fig.	BLC0B	BLC0B	BLC0B	BLC0B	BLC0B
Größe Size ² Tamaño	070	070	070	070	070
Länge Length	9,0	10,0	11,0	12,0	13,0
Longitud mm					
Bohrtiefe Drilling depth	5,0	4,0	3,0	2,0	1,0
Profundidad de perforación mm					

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
E-mail: info@meisingerusa.com

Balloon-Lift-Control

Balloon-Assisted
Sinus Floor Elevation

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Balloon-Lift-Control (BLC) is a very safe procedure used for the pre-augmentative stripping of the maxillary sinus mucosa. It is minimally invasive and extremely easy to pick up. Independently of the vertical height of the osseous sinus floor, the mucosa can be lifted by more than 10 mm. BLC consists of two osteotomes with the help of which a balloon expander is being inserted into the maxillary sinus after impression of the weakened floor. Subsequently the balloon is extended to the required size using radio-opaque liquid. The necessary volume equals the volume of augmentation required. This way, sinus floor augmentation becomes a procedure that can be carried out easily and quickly in all dental practices that specialise in implants without leading to the typical complications associated with conventional sinus lifts.

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The use of the initial bur combined with the respective shell prevents the special bur from slipping. If required, the allport bur or one of the abrasive burs can refresh and smooth the bone.

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Depending on the existing bone and the size of the implant to be inserted at a later date, now is the time to choose the guidance shell, the special bur and the osteotome with the larger (6.0 mm) or smaller (3.8 mm) diameter. Now the appropriate guidance shell is fastened onto the guidance tool, which is then held in front of the prepared alveolar crest.

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Subsequently and in order to thin the sinus floor you use the special bur to drill through the guidance shell into the alveolar crest until the thickness of the remaining bone measures 1 mm. The perfectly matching tools ensure that this pre-determined height of the remaining bone is maintained.

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The tip of the osteotome, with the help of the laser markings (mm scale), is aligned and set to the appropriate drill canal depth then fixed with the wider counter screw. The exit of the inner mandrin is, with help of the lower knurl-screw, adapted to the remaining bone-height accordingly. This is due to the help of the laser markings (mm scale) on the lower part of the mandrin with the subsequent fixation of the lower screw. With this one can precisely set the withdrawal point of the mandrin out of the osteotome. Subsequently the tip of the osteotome is positioned in the pre-formed canal the mandrel is then initiated.

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By slightly tapping on the mandrin grip, the compacta on the maxillary sinus side is fractured in a controlled fashion. Prior to this, the part of the mandrin that is now rising above the top of the shaft was accurately adjusted to the height of the remaining bone.

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Removal of the mandrin from the osteotome. Prior to inserting the balloon expander the air must be evacuated i.e. it must be pre-filled with fluid (see the detailed description of the balloon expander), then it is inserted into the guiding canal of the osteotome as far as it will go. In doing so, the blue extension at the end of the catheter shaft is pushed right down into the deepening of the osteotome grip.

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This is followed by a repeated (carried out at least 5 times) slow and controlled extension of the balloon with a maximum of approximately 3 ml of fluid (Ultravist® 240 or physiological, 0.9 % NaCl solution). This procedure ensures the stripping and lifting of the maxillary sinus mucosa. If required, a respective X-ray image can be generated with the extended balloon positioned in situ.

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At the end, the balloon is emptied again and removed with the osteotome. Then a spoon is used to insert the exact volume of the augmentation, which before served as the fluid to extend the balloon, step by step through the drilled hole into the newly established sinus area. By the application of the Osteotoms with the diameter of 3.8 mm, a direct implantation is possible.

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General instructions:

- All instruments are supplied non-sterile and must therefore be sterilised before use. Detailed instructions may be found in the MEISINGER Instructions and safety precautions leaflet for the use of rotary instruments in the medical area.

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